

Project Name:	RECOLOR: A Mobile Vision Assistive Tool for Color Vision Deficiency Using Adaptive Daltonization and Clustering Algorithms		
Written By:	CAPARAS, Mark; RULL Charles Justine; CEPEDA Johann Nikkolai; LACSA Stephen		
Written Date:	March 18 2026		
Test Case Description:	Functionality Test Case		
Testers:	Test Respondent's Full Name with Signature		
Test Conditions	Total Test Case:	233	Date:
	Pass:		Rating:
	Fail:		

Test Case Scenario ID	Name of the Module Function	Test Case Scenario	Action	Actual Input	Pass	Fail	Comments / Suggestions
AUTHENTICATION							
	<i>Sign Up</i>						
1. RC-AU	Sign Up – Open	Tapping Sign Up opens registration screen.	Registration fields displayed.	Tap Sign Up.			
2. RC-AU	Sign Up – Blank	All fields blank, tap submit.	Error: fields required.	Leave blank. Tap Submit.			
3. RC-AU	Sign Up – Invalid Email	Invalid email format entered.	Error: 'Invalid email.'	Enter 'abc@'. Tap Submit.			
4. RC-AU	Sign Up – No Password	Email entered, password blank.	Error: 'Password required.'	Enter email. Leave password blank.			
5. RC-AU	Sign Up – Short Pass	Password fewer than 6 chars.	Error: 'Password too short.'	Enter 5-char password.			

6. RC-AU	Sign Up – Valid	Valid credentials entered.	Account created.	Enter valid email + password.			
7. RC-AU	Sign Up – Duplicate	Already-registered email used.	Error: 'Email already in use.'	Enter existing email.			
	Sign In						
8. RC-AU	Sign In – Open	Tap Sign In from login.	Login fields displayed.	Tap Sign In.			
9. RC-AU	Sign In – Correct	Correct credentials entered.	Home screen displayed.	Enter correct email + pass.			
10. RC-AU	Sign In – Wrong Pass	Wrong password entered.	Error displayed. No crash.	Enter wrong password.			
11. RC-AU	Sign In – Unregistered	Unregistered email entered.	Error: 'Account not found.'	Enter unregistered email.			
12. RC-AU	Sign In – Blank Email	Email blank, tap sign in.	Error displayed.	Leave email blank.			
13. RC-AU	Sign In – Blank Pass	Password blank, tap sign in.	Error displayed.	Leave password blank.			
	Forgot Password						
14. RC-AU	Forgot Pass – Open	Tap Forgot Password.	Reset screen opens.	Tap Forgot Password.			
15. RC-AU	Forgot Pass – Valid	Registered email entered.	'Reset email sent.'	Enter registered email.			

16. RC-AU	Forgot Pass – Invalid	Unregistered email entered.	Error or info displayed.	Enter unregistered email.			
	<i>Sign Out & Session</i>						
17. RC-AU	Sign Out	Tap Sign Out.	Login screen shown.	Tap Sign Out.			
18. RC-AU	Sign Out – Back	Back button after sign out.	Remains on login.	Press back after sign out.			
19. RC-AU	Session – Restart App	Close and reopen app.	Still signed in.	Close app. Reopen.			
20. RC-AU	Session – Restart Device	Restart device, open app.	Still signed in.	Restart device. Open app.			
HOME SCREEN AND NAVIGATION							
	<i>Home Screen</i>						
21. RC-NAV	Home Screen	All main features visible.	Camera, Identifier, Ishihara, Career, Settings shown.	Observe home after login.			
22. RC-NAV	Icons/Labels	All icons labeled or identifiable.	No unlabeled icons.	Observe all icons.			
	<i>Screen Navigation</i>						
23. RC-NAV	Open Camera	Tap Camera feature.	Camera screen opens.	Tap Camera.			
24. RC-NAV	Open Identifier	Tap Color Identifier.	Identifier screen opens.	Tap Identifier.			
25. RC-NAV	Open Ishihara	Tap Ishihara Screening.	Ishihara module opens.	Tap Ishihara.			

26. RC-NAV	Open Career	Tap Career Awareness.	Career screen opens.	Tap Career.			
27. RC-NAV	Open Settings	Tap Settings.	Settings screen opens.	Tap Settings.			
	Back Button						
28. RC-NAV	Back – Camera	Back from Camera.	Home shown.	Press back.			
29. RC-NAV	Back – Identifier	Back from Identifier.	Home shown.	Press back.			
30. RC-NAV	Back – Ishihara	Back from Ishihara.	Home shown.	Press back.			
31. RC-NAV	Back – Career	Back from Career.	Home shown.	Press back.			
32. RC-NAV	Back – Home	Back from Home.	App exits or confirmation.	Press back.			
CAMERA MODULE							
	Permissions						
33. RC-CAM	Permission Dialog	First time camera open shows permission.	Permission dialog appears.	Open camera (fresh install).			
34. RC-CAM	Permission Granted	Grant permission.	Live feed displays.	Tap Allow.			
35. RC-CAM	Permission Denied	Deny permission.	Message shown. No crash.	Tap Deny.			
36. RC-CAM	Permission Permanent	Deny permanently.	Directs to settings.	Tap Don't Ask Again.			
	Rear Camera Feed						
37. RC-CAM	Rear Live Feed	Rear camera shows smooth feed.	No black screen.	Open camera.			

38. RC-CAM	Rear No Freeze (2 min)	Feed runs 2+ minutes.	No freeze, no ANR.	Keep open 2 minutes.			
39. RC-CAM	Rear Extended (10 min)	Feed runs 10+ minutes.	No freeze, no memory leak.	Keep open 10 minutes.			
40. RC-CAM	Background/Foreground	Background and return.	Feed resumes.	Background. Return.			
41. RC-CAM	Phone Call Resume	Receive call and return.	Feed resumes.	Receive call. Return.			
42. RC-CAM	Lock/Unlock	Lock and unlock phone.	Feed resumes.	Lock. Unlock. Return.			
43. RC-CAM	Low Light	Camera in low-light.	Feed displays (dark, no crash).	Cover lens partially.			
	Front Camera						
44. RC-CAM	Front Camera Feed	Front camera shows live feed.	Selfie view displayed.	Switch to front camera.			
45. RC-CAM	Switch Rear→Front	Tap camera switch button.	Front camera activates.	Tap switch button.			
46. RC-CAM	Switch Front→Rear	Tap camera switch again.	Rear camera activates.	Tap switch button.			
47. RC-CAM	Correction on Front	Color correction works on front camera.	Correction visible on selfie feed.	Enable correction. Switch to front.			

48. RC-CAM	Identifier on Front	Color Identifier works on front camera.	Label displays for front feed.	Switch to front. Open Identifier.			
49. RC-CAM	Simulation on Front	CVD Simulation works on front camera.	Simulation visible on front feed.	Enable simulation. Switch to front.			
50. RC-CAM	Rapid Switch (5x)	Switch front/rear 5 times rapidly.	No crash. Feed recovers.	Toggle switch 5 times.			
	<i>Shutter / Capture</i>						
51. RC-CAM	Shutter Button	Tap shutter to capture.	Image captured. Preview shown.	Tap shutter button.			
52. RC-CAM	Capture + Correction ON	Capture with correction active.	Saved image shows corrected colors.	Enable correction. Tap shutter.			
53. RC-CAM	Capture + Correction OFF	Capture with correction off.	Saved image shows original colors.	Disable correction. Tap shutter.			
54. RC-CAM	Capture Preview	Preview/confirmation shown after capture.	User can accept/discard or auto-save.	Tap shutter. Observe.			
55. RC-CAM	Save to Device	Captured image saves to storage.	Image appears in device gallery.	Capture. Check device gallery.			
	<i>Freeze / Pause</i>						
56. RC-CAM	Freeze Button	Tap freeze/pause button.	Camera feed pauses on current frame.	Tap freeze button.			

66. RC-CAM	Switch: Identifier→Correction	Switch from Identifier to Correction.	Feed continues. No blackout.	Switch mode.			
67. RC-CAM	Switch: Correction→Simulation	Switch from Correction to Simulation.	Feed continues. No blackout.	Switch mode.			
68. RC-CAM	Switch: Simulation→Identifier	Switch from Simulation to Identifier.	Feed continues. No blackout.	Switch mode.			
69. RC-CAM	Rapid Mode Switch	Switch all 3 modes rapidly.	No crash. No permanent blackout.	Toggle modes 5+ times.			
70. RC-CAM	Blackout Recovery	If blackout occurs, swap front/rear recovers.	Feed restored after camera swap.	If blackout: tap camera switch.			
71. RC-CAM	Mode Switch 10+ Times	Switch modes 10+ times consecutively.	No crash, no memory leak.	Cycle modes 10 times.			
COLOR CORRECTION (DALTONIZATION)							
	<i>Correction Toggle</i>						
72. RC-CC	Enable Correction	Enable correction.	Visible color change.	Toggle ON.			
73. RC-CC	Disable Correction	Disable correction.	Normal colors.	Toggle OFF.			
74. RC-CC	Rapid Toggle (5x)	Toggle 5 times rapidly.	No crash/freeze.	Toggle 5x quickly.			

75. RC-CC	Correction While Moving	Move camera with correction ON.	Continuous correction.	Move camera around.			
	<i>CVD Type Selector</i>						
76. RC-CC	Select Protan	Select Protanomaly.	Protan correction.	Select Protanomaly.			
77. RC-CC	Select Deutan	Select Deuteranomaly.	Deutan correction.	Select Deuteranomaly.			
78. RC-CC	Switch Pro→Deu	Switch while ON.	Visible change.	Switch type.			
79. RC-CC	Switch Deu→Pro	Switch while ON.	Visible change.	Switch type.			
80. RC-CC	Type Persists	Leave and return to camera.	Same type.	Select. Leave. Return.			
	<i>Strength Slider</i>						
81. RC-CC	Slider Visible (ON)	Slider shown when ON.	Slider displayed.	Enable correction.			
82. RC-CC	Slider Hidden (OFF)	Slider hidden when OFF.	Not interactive.	Disable correction.			
83. RC-CC	Slider 0.0	Drag to minimum.	No correction.	Drag to 0.0.			
84. RC-CC	Slider 0.25	Drag to 0.25.	Mild correction.	Drag to 0.25.			
85. RC-CC	Slider 0.5	Drag to midpoint.	Moderate correction.	Drag to 0.5.			
86. RC-CC	Slider 0.75	Drag to 0.75.	Strong correction.	Drag to 0.75.			
87. RC-CC	Slider 1.0	Drag to maximum.	Full correction.	Drag to 1.0.			
88. RC-CC	Proportional	Smooth change 0→1.	No jumps.	Drag slowly.			

89. RC-CC	Rapid Drag	Drag rapidly.	No crash/tearing.	Drag back/forth.			
90. RC-CC	Persist (Screen)	Leave and return.	Same position.	Set 0.5. Leave. Return.			
91. RC-CC	Persist (BG)	Background and return.	Same position.	Set 0.5. BG. Return.			
92. RC-CC	Type Switch at 0.5	Switch type at 0.5.	New type, same strength.	Set 0.5. Switch type.			
	Correction Quality						
93. RC-CC	Red/Green Deutan	Deutan on red/green scene.	More distinguishable.	Point at red/green.			
94. RC-CC	Red/Green Protan	Protan on red/green scene.	More distinguishable.	Point at red/green.			
95. RC-CC	No Bleeding	No color bleeding.	Clean boundaries.	Observe boundaries.			
96. RC-CC	No Flicker	No flickering.	Stable output.	Watch 30 seconds.			
97. RC-CC	No Banding	No banding.	Smooth gradients.	Point at gradient.			
98. RC-CC	Extended (5 min)	5+ minutes use.	No degradation.	Keep ON 5 min.			
CVD SIMULATION MODE							
99. RC-SIM	Enable Simulation	Enable simulation.	Colors shift.	Toggle ON.			
100. RC-SIM	Disable Simulation	Disable simulation.	Normal colors.	Toggle OFF.			
101. RC-SIM	Protan Simulation	Protan reduces red.	Reds olive/brown.	Select Protan.			
102. RC-SIM	Deutan Simulation	Deutan reduces green.	Greens brown/yellow.	Select Deutan.			

103. RC-SIM	Sim $\neq$ Correction	Different from correction.	Two distinct effects.	Compare modes.			
104. RC-SIM	Toggle Sim/Corr	Switch between modes.	Clear difference.	Switch modes.			
105. RC-SIM	Multi-Device	Works on all 3 devices.	Consistent.	Test all devices.			
COLOR IDENTIFIER							
	<i>10-Color Identification</i>						
106. RC-CI	Red	Point at red.	Label: 'Red.'	Point at red.			
107. RC-CI	Green	Point at green.	Label: 'Green.'	Point at green.			
108. RC-CI	Blue	Point at blue.	Label: 'Blue.'	Point at blue.			
109. RC-CI	Yellow	Point at yellow.	Label: 'Yellow.'	Point at yellow.			
110. RC-CI	Orange	Point at orange.	Label: 'Orange.'	Point at orange.			
111. RC-CI	Pink	Point at pink.	Label: 'Pink.'	Point at pink.			
112. RC-CI	Brown	Point at brown.	Label: 'Brown.'	Point at brown.			
113. RC-CI	Violet	Point at violet.	Label: 'Violet.'	Point at violet.			
114. RC-CI	Cyan	Point at cyan.	Label: 'Cyan.'	Point at cyan.			
115. RC-CI	Neutral	Point at white/gray/black.	Label: 'Neutral.'	Point at neutral.			
	<i>Responsiveness</i>						
116. RC-CI	Red→Blue	Move red to blue.	Updates ~1s.	Move camera.			

117. RC-CI	Green→Yellow	Move green to yellow.	Updates ~1s.	Move camera.			
118. RC-CI	Rapid Movement	Rapid color changes.	No crash.	Move rapidly.			
119. RC-CI	Steady 30s	Hold on one color.	Stable label.	Hold 30 seconds.			
	<i>Edge Cases</i>						
120. RC-CI	Multi-Color	Multi-colored scene.	Dominant color.	Point at mixed.			
121. RC-CI	Very Dark	Near-black scene.	No crash.	Point at dark.			
122. RC-CI	Very Bright	Near-white scene.	No crash.	Point at bright.			
123. RC-CI	Lens Covered	Cover lens.	No crash.	Cover lens.			
	<i>TFLite Model</i>						
124. RC-CI	TFLite Load	Model loads.	No crash.	Open Identifier.			
125. RC-CI	Multi-Device	Works on 3 devices.	Labels on all.	Test all devices.			
ISHIHARA SCREENING MODULE							
	<i>Launch and Setup</i>						
126. RC-ISH	Open Module	Open Ishihara.	Module opens.	Tap Ishihara.			
127. RC-ISH	Brightness Check	Brightness < 80%.	Prompt shown.	Set 50%. Open.			
128. RC-ISH	Brightness OK	Brightness ≥ 80%.	Plate 1 shown.	Set 100%. Open.			
	<i>Plate Display (1–25)</i>						
129. RC-ISH	Plate 1 (Demo)	Plate renders correctly.	Correct colors. No errors.	View Plate 1.			

130. RC-ISH	Plate 2 (Transform)	Plate renders correctly.	Correct colors. No errors.	View Plate 2.			
131. RC-ISH	Plate 3 (Transform)	Plate renders correctly.	Correct colors. No errors.	View Plate 3.			
132. RC-ISH	Plate 4 (Transform)	Plate renders correctly.	Correct colors. No errors.	View Plate 4.			
133. RC-ISH	Plate 5 (Transform)	Plate renders correctly.	Correct colors. No errors.	View Plate 5.			
134. RC-ISH	Plate 6 (Transform)	Plate renders correctly.	Correct colors. No errors.	View Plate 6.			
135. RC-ISH	Plate 7 (Transform)	Plate renders correctly.	Correct colors. No errors.	View Plate 7.			
136. RC-ISH	Plate 8 (Transform)	Plate renders correctly.	Correct colors. No errors.	View Plate 8.			
137. RC-ISH	Plate 9 (Transform)	Plate renders correctly.	Correct colors. No errors.	View Plate 9.			
138. RC-ISH	Plate 10 (Vanish)	Plate renders correctly.	Correct colors. No errors.	View Plate 10.			
139. RC-ISH	Plate 11 (Vanish)	Plate renders correctly.	Correct colors. No errors.	View Plate 11.			
140. RC-ISH	Plate 12 (Vanish)	Plate renders correctly.	Correct colors. No errors.	View Plate 12.			

141. RC-ISH	Plate 13 (Vanish)	Plate renders correctly.	Correct colors. No errors.	View Plate 13.			
142. RC-ISH	Plate 14 (Vanish)	Plate renders correctly.	Correct colors. No errors.	View Plate 14.			
143. RC-ISH	Plate 15 (Vanish)	Plate renders correctly.	Correct colors. No errors.	View Plate 15.			
144. RC-ISH	Plate 16 (Vanish)	Plate renders correctly.	Correct colors. No errors.	View Plate 16.			
145. RC-ISH	Plate 17 (Vanish)	Plate renders correctly.	Correct colors. No errors.	View Plate 17.			
146. RC-ISH	Plate 18 (Hidden)	Plate renders correctly.	Correct colors. No errors.	View Plate 18.			
147. RC-ISH	Plate 19 (Hidden)	Plate renders correctly.	Correct colors. No errors.	View Plate 19.			
148. RC-ISH	Plate 20 (Hidden)	Plate renders correctly.	Correct colors. No errors.	View Plate 20.			
149. RC-ISH	Plate 21 (Hidden)	Plate renders correctly.	Correct colors. No errors.	View Plate 21.			
150. RC-ISH	Plate 22 (Classify)	Plate renders correctly.	Correct colors. No errors.	View Plate 22.			
151. RC-ISH	Plate 23 (Classify)	Plate renders correctly.	Correct colors. No errors.	View Plate 23.			

152. RC-ISH	Plate 24 (Classify)	Plate renders correctly.	Correct colors. No errors.	View Plate 24.			
153. RC-ISH	Plate 25 (Classify)	Plate renders correctly.	Correct colors. No errors.	View Plate 25.			
	<i>User Input</i>						
154. RC-ISH	Input Response	Can input per plate.	Field/buttons available.	Enter response.			
155. RC-ISH	Advance Plate	Can advance to next.	Next plate loads.	Tap Next.			
156. RC-ISH	Plate Order	Sequential order.	1→21, then 22–25.	Navigate all.			
	<i>Scoring — Stage 1</i>						
157. RC-ISH	21/21 → Normal	All correct.	Normal.	Enter all correct.			
158. RC-ISH	17/21 → Normal	Boundary.	Normal.	Enter 17 correct.			
159. RC-ISH	16/21 → Indeterminate	Boundary.	Indeterminate.	Enter 16 correct.			
160. RC-ISH	14/21 → Indeterminate	14 correct.	Indeterminate.	Enter 14 correct.			
161. RC-ISH	13/21 → CVD	Boundary.	CVD. Stage 2.	Enter 13 correct.			
162. RC-ISH	5/21 → Severe	5 correct.	Severe.	Enter 5 correct.			
163. RC-ISH	0/21	Zero correct.	No crash.	Enter 0 correct.			
	<i>Scoring — Stage 2</i>						
164. RC-ISH	Protan Pattern	Protan answers.	Protanomaly.	Enter 6,2,5,6.			
165. RC-ISH	Deutan Pattern	Deutan answers.	Deuteranomaly.	Enter 2,4,3,9.			
166. RC-ISH	Mixed Pattern	Mixed 2+2.	Indeterminate.	Enter mixed.			

167. RC-ISH	Normal Pattern	Normal answers.	Handles gracefully.	Enter 26,42,35,96.			
	<i>Pre/Post Flow</i>						
168. RC-ISH	Pre-Test Score	Complete pre-test.	Score saved.	Complete without correction.			
169. RC-ISH	Enable Post-Corr	Enable correction after pre.	Correction activates.	Enable correction.			
170. RC-ISH	Post-Test Score	Complete post-test.	Score saved.	Complete with correction.			
171. RC-ISH	Score Comparison	Both scores shown.	Improvement calculated.	Observe results.			
172. RC-ISH	Enforce Order	Cannot skip pre to post.	Post unavailable.	Try post without pre.			
HAPTIC FEEDBACK (VIBRATION)							
173. RC-HF	Haptic ON	Setting ON.	Reflects ON.	Check setting.			
174. RC-HF	Red→Green	Color change.	Vibration felt.	Point red→green.			
175. RC-HF	Blue→Yellow	Color change.	Vibration felt.	Point blue→yellow.			
176. RC-HF	Steady	One color held.	No vibration.	Hold steady.			
177. RC-HF	Toggle OFF	Disable haptic.	OFF reflected.	Toggle OFF.			
178. RC-HF	OFF: No Vibration	Haptic OFF, change colors.	No vibration.	Move colors.			
179. RC-HF	Toggle Back ON	Re-enable.	Resumes.	Toggle ON.			
180. RC-HF	Silent Mode	Phone silent.	Still vibrates.	Set silent. Test.			

181. RC-HF	Multi-Device	All 3 devices.	Consistent.	Test all.			
AUDIO / VOICE FEEDBACK (TTS)							
182. RC-AF	Audio ON	Setting ON.	Reflects ON.	Check setting.			
183. RC-AF	Speak Red	Point at red.	'Red' spoken.	Point at red.			
184. RC-AF	Speak Blue	Point at blue.	'Blue' spoken.	Point at blue.			
185. RC-AF	Update	Move colors.	New label ~1–2s.	Move red→green.			
186. RC-AF	Rapid	Rapid changes.	No overlap.	Move rapidly.			
187. RC-AF	Toggle OFF	Disable audio.	OFF reflected.	Toggle OFF.			
188. RC-AF	OFF: No Audio	Audio OFF.	No audio.	Point at colors.			
189. RC-AF	Toggle Back ON	Re-enable.	Resumes.	Toggle ON.			
190. RC-AF	Volume	Device volume.	Quieter.	Lower volume.			
191. RC-AF	Headphones	Headphones connected.	Audio in headphones.	Connect headphones.			
192. RC-AF	Multi-Device	All 3 devices.	Consistent.	Test all.			
CAREER AWARENESS MODULE							
193. RC-CA	Open Screen	Open from home.	Screen loads.	Tap Career.			
194. RC-CA	Text Displayed	Text fully shown.	No truncation.	Observe text.			
195. RC-CA	Images Load	Graphics load.	No broken images.	Observe graphics.			
196. RC-CA	Scrollable	Content scrollable.	Can scroll all.	Scroll down.			
197. RC-CA	Professions Listed	CVD professions shown.	At least 3.	Read content.			

198. RC-CA	Accurate Info	Factually accurate.	No errors.	Verify content.			
199. RC-CA	Back Button	Back returns.	Previous screen.	Press back.			
200. RC-CA	Multi-Device	All 3 devices.	No issues.	Test all.			
SETTINGS AND USER PROFILE							
	CVD Profile						
201. RC-SET	Select CVD Type	Select Protan/Deutan/None.	Saved.	Select type.			
202. RC-SET	Type in Correction	Reflected in correction.	Uses selected.	Open camera.			
203. RC-SET	Type Persists	Persists after restart.	Same type.	Restart. Check.			
	Toggles						
204. RC-SET	Haptic Toggle	Present and works.	ON/OFF works.	Toggle haptic.			
205. RC-SET	Audio Toggle	Present and works.	ON/OFF works.	Toggle audio.			
206. RC-SET	Independent	One doesn't affect other.	Independent.	Toggle one. Check other.			
	Data Persistence						
207. RC-SET	Profile Display	Profile info shown.	Name/email shown.	Open profile.			
208. RC-SET	Profile Persists	Persists after restart.	Not lost.	Restart. Check.			
209. RC-SET	Settings Persist	Settings persist.	Maintained.	Restart. Check.			
FIREBASE AND DATA STORAGE							
	Firestore Logging						

210. RC-FB	Pre-Test Logged	Data in Firebase.	All fields populated.	Complete pre. Check Firebase.			
211. RC-FB	Post-Test Logged	Post data logged.	correction_active=true.	Complete post. Check.			
212. RC-FB	Per-Plate Data	Individual plates stored.	plate_number, response, etc.	Check fields.			
213. RC-FB	CVD Type Logged	Classification logged.	Type recorded.	Check field.			
	Data Integrity						
214. RC-FB	Score Match	Firebase = app display.	No discrepancy.	Compare.			
215. RC-FB	Timestamp	Accurate timestamp.	Within 1 minute.	Check time.			
216. RC-FB	No Duplicates	One record per test.	No duplicates.	Check count.			
	Security						
217. RC-FB	No Sensitive Logs	logcat clean.	No CVD/personal data.	Run logcat.			
218. RC-FB	HTTPS Only	All traffic HTTPS.	No plaintext.	Monitor network.			
219. RC-FB	Auth Required	Unauth blocked.	Permission denied.	Test direct API.			
	Offline						
220. RC-FB	Offline Test	Test without WiFi.	No crash.	Disable WiFi. Test.			
221. RC-FB	Offline Sync	Reconnect synes.	Data in Firebase.	Re-enable WiFi.			
GALLERY AND IMAGE MANAGEMENT							
	Gallery Screen						

222. RC-GAL	Open Gallery	Open from home/menu.	Gallery screen loads.	Tap Gallery.			
223. RC-GAL	Display Images	Shows saved images.	Thumbnails visible.	Observe gallery.			
224. RC-GAL	Full View	Tap image for full size.	Full image displayed.	Tap thumbnail.			
225. RC-GAL	Before/After	Comparison view available.	Original + corrected shown.	Tap compare option.			
226. RC-GAL	Empty State	Gallery with no images.	'No images yet' message.	Open with no saves.			
	Upload & Process						
227. RC-GAL	Upload Image	Upload from device library.	Image picker opens. Image loads.	Tap upload/import.			
228. RC-GAL	Correction on Upload	Apply correction to uploaded image.	Corrected version displayed.	Upload. Apply correction.			
229. RC-GAL	Simulation on Upload	Apply simulation to uploaded image.	Simulated version displayed.	Upload. Apply simulation.			
230. RC-GAL	Identifier on Upload	Apply identifier to uploaded image.	Color label shown.	Upload. Apply identifier.			
	Save & Share						
231. RC-GAL	Save Corrected	Save corrected image to device.	Image in device gallery.	Tap save.			

232. RC-GAL	Storage Permission	Permission requested for save/upload.	Permission dialog shown.	First save/upload.			
233. RC-GAL	Permission Denied	Deny storage permission.	Message shown. No crash.	Tap Deny.			

Total Test Criteria: 233

Passed: _____ / 233

Failed: _____ / 233

Pass Rate: _____ %

Failed Criteria (list):

Comments / Recommendations:

Test Respondent's Signature Over Printed Name / Lagda ng Tagasuri

Date / Petsa: _____

Project Name:	RECOLOR: A Mobile Vision Assistive Tool for Color Vision Deficiency Using Adaptive Daltonization and Clustering Algorithms			
Written By:				
Written Date:				
Test Case Description:	Compatibility Test Case			
Testers:	Test Respondent's Full Name with Signature			
Test Conditions	Total Test Case:	19	Date:	
	Pass:		Rating:	
	Fail:			

Test Case Scenario ID	Name of the Module Function	Test Case Scenario	Action	Actual Input	Pass	Fail	Comments / Suggestions
	Compatibility						
1. RC-GAL	Multi-Device	Gallery works on all 3 devices.	No issues.	Test all devices.			
COMPATIBILITY TESTING							
	Device 1 (Low-end): _____						
2. RC-COM	Device 1	Installs.	Installs	Test on Device 1 (Low-end).			
3. RC-COM	Device 1	Launches.	Launches	Test on Device 1 (Low-end).			
4. RC-COM	Device 1	Camera works.	Camera works	Test on Device 1 (Low-end).			
5. RC-COM	Device 1	Correction works.	Correction works	Test on Device 1 (Low-end).			

6. RC-COM	Device 1	Ishihara displays.	Ishihara displays	Test on Device 1 (Low-end).			
	Device 2 (Mid-range): _____						
7. RC-COM	Device 2	Installs.	Installs	Test on Device 2 (Mid-range).			
8. RC-COM	Device 2	Launches.	Launches	Test on Device 2 (Mid-range).			
9. RC-COM	Device 2	Camera works.	Camera works	Test on Device 2 (Mid-range).			
10. RC-COM	Device 2	Correction works.	Correction works	Test on Device 2 (Mid-range).			
11. RC-COM	Device 2	Ishihara displays.	Ishihara displays	Test on Device 2 (Mid-range).			
	Device 3 (High-end): _____						
12. RC-COM	Device 3	Installs.	Installs	Test on Device 3 (High-end).			
13. RC-COM	Device 3	Launches.	Launches	Test on Device 3 (High-end).			
14. RC-COM	Device 3	Camera works.	Camera works	Test on Device 3 (High-end).			
15. RC-COM	Device 3	Correction works.	Correction works	Test on Device 3 (High-end).			

16. RC-COM	Device 3	Ishihara displays.	Ishihara displays	Test on Device 3 (High-end).			
	Display						
17. RC-COM	Smallest Screen	UI correct on smallest.	No overlaps.	Inspect smallest.			
18. RC-COM	Largest Screen	UI correct on largest.	Scales properly.	Inspect largest.			
19. RC-COM	Touch Targets	Buttons usable all sizes.	Tappable.	Test all.			

Total Test Criteria: 19

Passed: _____ / 19

Failed: _____ / 19

Pass Rate: _____ %

Failed Criteria (list):

Comments / Recommendations:

Test Respondent's Signature Over Printed Name / Lagda ng Tagasuri

Date / Petsa: _____